



Indication

1. High-flow irrigation in arthroscopic approach procedures.
2. Irrigation in urological surgeries, with effective control of bleeding.
3. Irrigation of cavities and/or perfusion of organs in surgery.

Description

1. Crystal PVC equipment with 2 or 4 extension routes, with diameters of 1/4".
2. Each route has a piercing tip, cap and flow clamp.
3. Crystal PVC connector at the end, with the option of a luer lock connector.
4. Latex-free product.
5. Sterilized in ETO.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
81003	ARTROZAMM 2 WAY	S, 1/4" x 1,6mm x 2,0m	10216350113	25
81006	ARTROZAMM 4 WAY	S, 1/4" x 1,6mm x 2,0m	10216350113	15
81007	ARTROZAMM 4 WAY	S, 1/4" x 1,6mm x 2,0m	10216350113	15

S Sterile

NS Non-sterile

V With Valve

Products subject to change without prior notice. Please contact us for availability.