



## Indication

1. Vacuum aspiration of airway secretions.

## Description

1. System suitable for aspirating secretions.
2. Available in 500 ml, 1000 ml and 2000 ml capacities.
3. Latex-free, minimising the risk of allergies.
4. Sterilised in ETO or available in a non-sterile version.
5. Integrated safety valve to prevent overflow.
6. The bottle complies with ABNT NBR 12188:2016 standard (use between 11.81 inHg and 18.11 inHg) for vacuum pressure.

| CODE | PRODUCT                               | DESCRIPTION                                         | ANVISA      | BOX |
|------|---------------------------------------|-----------------------------------------------------|-------------|-----|
| 1425 | KIT AVAZAMM 2000ml AR PLUS (LIPO KIT) | NS, V, 1/4" x 1,6mm x 3,0m, 1/4" x 1,6mm x 3,0m     | 10216350070 | 10  |
| 1440 | KIT AVAZAMM 1000ML AR PLUS            | S, V, 3/16" x 1,0mm x 1,20m, 3/16" x 1,0mm x 0,50m  | 10216350091 | 20  |
| 1444 | KIT AVAZAMM 1000ML AR PLUS            | NS, V, 3/16" x 1,0mm x 1,20m, 3/16" x 1,0mm x 0,50m | 10216350070 | 25  |
| 1447 | KIT AVAZAMM 1000ML 2 VIAS AR PLUS     | S, V, 1/4" x 1,6mm x 1,50m, 3/16" x 1,0mm x 1,50m   | 10216350091 | 15  |
| 1485 | KIT AVAZAMM 1000ML AR PLUS            | S, V, 3/16" x 1,0mm x 2,50m, 3/16" x 1,0mm x 2,50m  | 10216350091 | 15  |
| 1492 | KIT AVAZAMM 500ML AR PLUS             | NS, V, 3/16" x 1,0mm x 1,20m, 3/16" x 1,0mm x 0,50m | 10216350070 | 30  |
| 1495 | KIT AVAZAMM 500ML AR PLUS             | S, V, 3/16" x 1,0mm x 1,20m, 3/16" x 1,0mm x 0,50m  | 10216350091 | 25  |
| 1456 | KIT AVAZAMM 1000ML AR PLUS            | S, V, 3/16" x 1,0mm x 1,90m, 3/16" x 1,0mm x 0,50m  | 10216350091 | 20  |

*S* Sterile

*NS* Non-sterile

*V* With Valve

Products subject to change without prior notice. Please contact us for availability.