

PREMERE

Pressure Bag

The solution for fluid infusion.

- Adjustable diameter bag for pressurizing flexible bottles or saline solution bags.
- Stronger outer plastic film and thinner inner film allow for faster and more consistent inflation. The bag wraps around the serum bottle at 360°, allowing for more uniform pressurization.
- Chrome pressure gauge with 60mm diameter and easy-to-view customized scale.
- Velcro fasteners, positioned so as not to obstruct the label view, facilitating the installation of the serum bottle in the bag.
- Large inflation bulb and caliber tube, for faster inflation and deflation.
- Durable and easy to clean.
- Available in 500ml and 1000ml.



Code	Description	Anvisa	Box
7510	Bolsa Pressórica PREMERE 500ml	10216350119	1 / 6 / 16 / 48
7511	Bolsa Pressórica PREMERE 500ml	10216350119	1 / 4 / 8 / 24

ZAMMI INSTRUMENTAL EIRELI

R. Bernardo Vasconcelos, 992 Parque Santa Lúcia, Duque de Caxias - RJ CEP: 25251-300

Tel: +55 21 2677-1500 | +55 11 5574-5767 | faleconosco@zammi.com.br

www.zammi.com.br