



Indication

1. Collection of tracheal secretion and/or bronchoalveolar lavage by direct aspiration or bronchoscopy.

Description

1. Crystal PVC collection bottle, graduated every 10 ml.
2. 120 ml capacity.
3. Comes with adhesive label to identify the patient and the sample collected.
4. Latex-free.
5. Sterilized in ETO.

| CODE | PRODUCT | DESCRIPTION | ANVISA | BOX |
|------|------------|-----------------------------|-------------|-----|
| 1451 | BRONCOZAMM | S, Broncozamm 120ml e 0,15m | 10216350096 | 40 |

S Sterile

NS Non-sterile

V With Valve

Products subject to change without prior notice. Please contact us for availability.