



## Indication

1. Aspiration in large cavities, orthopedic surgeries and organ removal.
2. Conventional suction in the operating room.

## Description

1. Swivel or flexible suction nozzle with flow regulator.
2. 3 meter PVC extension with flexible connector.
3. Non-toxic and thermogenic material.
4. Latex-free, ensuring greater safety.
5. Sterilized in ETO or available in a non-sterile version

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
1230	DRENOSUCZAMM EXTRA 20/8	S, 20/8 x 3,0m   Ponteira rígida	10216350054	30
1222	DRENOSUCZAMM LONG SOFT 18/7	S, 18/7 x 3,0m   Ponteira flexível	10216350054	30

*S* Sterile

*NS* Non-sterile

*V* With Valve

Products subject to change without prior notice. Please contact us for availability.