




Indication

1. Drainage of the intrathoracic cavities (thoracic, pleural and mediastinal).

Description

1. Radiopaque chest drain.
2. Latex-free.
3. Sterilized in ETO.
4. Made of PVC in lengths of 0.40 m, 0.45 m and 0.50 m, available in gauges from 10 to 40 French.
5. It has a radiopaque filament throughout, an atraumatic tip and is easily controlled, allowing the drain to be clamped.
Includes optional connector for connection to the drainage system.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
2201	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 10	S, 10 Frenchs, 0,40m	10216350112	1/600
2205	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 12	S, 12 Frenchs, 0,40m	10216350112	1/600
2210	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 14	S, 14 Frenchs, 0,40m	10216350112	1/600
2216	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 16	S, 16 Frenchs, 0,40m	10216350112	1/600
2220	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 18	S, 18 Frenchs, 0,40m	10216350112	1/600
2225	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 20	S, 20 Frenchs, 0,40m	10216350112	1/300
2230	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 22	S, 22 Frenchs, 0,40m	10216350112	1/300
2235	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 24	S, 24 Frenchs, 0,45m	10216350112	1/300
2240	DTUBULAR CHEST DRAIN - DRENOZAMM LT FS 26	S, 26 Frenchs, 0,45m	10216350112	1/300
2245	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 28	S, 28 Frenchs, 0,45m	10216350112	1/300
2250	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 30	S, 30 Frenchs, 0,45m	10216350112	1/300
2255	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 32	S, 32 Frenchs, 0,45m	10216350112	1/300
2260	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 34	S, 34 Frenchs, 0,45m	10216350112	1/300
2265	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 36	S, 36 Frenchs, 0,45m	10216350112	1/300
2270	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 38	S, 38 Frenchs, 0,45m	10216350112	1/300
2275	TUBULAR CHEST DRAIN - DRENOZAMM LT FS 40	S, 40 Frenchs, 0,45m	10216350112	1/300

S Sterile

NS Non-sterile

V With Valve

Products subject to change without prior notice. Please contact us for availability.