



Indication

- 1. CO2 insufflator for videolaparoscopy.

Description

- 1. Non-toxic crystal PVC extension, 1/4" and 3m
- 2. length.
- 3. Filter for bacterial and viral barrier.
- 4. Latex-free.
- 5. Sterilized in ETO.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
8721	EXTENSÃO IVL	S, 1,4" x 1,6mm x 3,0m, smooth connector	10216350093	25

S Sterile
NS Non-sterile
V With Valve
Products subject to change without prior notice. Please contact us for availability.