



Indication

1. Simultaneous irrigation and aspiration in surgical procedures.

Description

1. 1/4" x 3m non-toxic PVC double extension, blue for transparency and crystal for suction.
2. 1/4" flexible PVC connector.
3. Flat choke for flow control.
4. Latex-free.
5. Sterilized in ETO.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
8200	IRRIGAZAMM DUO	S, 1/4" x 3,0m, with flexible connector, 2-way	10216359014	20

S Sterile

NS Non-sterile

V With Valve

Products subject to change without prior notice. Please contact us for availability.