



Indication

1. Procedure for aspirating secretions and body fluids during surgical procedures in a hospital environment.
2. Reduction in the volume of hospital waste and the amount of plastic discarded, as well as optimization of storage space.

Description

1. Reservoir available in 1,000ml, 2,000ml and 3,000ml capacities.
2. Anti-reflux valve with bacteriological filter, protecting the vacuum network and preventing overflow during transportation.
3. Latex-free, ensuring safety for sensitive users.
4. Non-sterile
5. Easy handling and safe disposal, minimizing the risk of contamination.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
1567	COLLECTOR SUPPORT MULTIZAMM ECO 2000ml	COLLECTOR. SUPPORT MULTIZAMM ECO 2000	10216350094	1
1592	COLLECTOR SUPPORT MULTIZAMM ECO 2000ml	NS, SV, GREEN LID	10216350094	1
1596	MULTIZAMM ECO 1000ml GREEN w/ valve	NS, V, GREEN LID	10216350094	25
1599	MULTIZAMM ECO 1000ML WHITE BOTTLE	Rigid bottle with wall bracket	10216350094	1
1507	MULTIZAMM ECO 2000ml	NS, SV, WHITE LID	10216350094	25
1508	MULTIZAMM ECO 2000ml GREEN w/ valve	NS, V, GREEN LID	10216350094	25
1562	MULTIZAMM ECO 3000ml GREEN W/ valve	NS, V, GREEN LID	10216350094	25
1563	MULTIZAMM ECO 3000ml WHITE	NS, V, WHITE LID	10216350094	25

S Sterile

NS Non-sterile

V With Valve

Products subject to change without prior notice. Please contact us for availability.