



## Indication

1. Procedure for aspirating secretions and body fluids during surgical procedures in a hospital environment.
2. Collection of secretions and purging through aspiration.

## Description

1. Rigid crystal PVC bottle, graduated 100/100ml, with a capacity of 2,000ml, 3,500ml and 5,000ml.
2. Quick screw cap and safety valve that stops aspiration when the bottle is full.
3. Anti-reflux valve with bacteriological filter, guaranteeing safety during use.
4. Two 2.5m crystal PVC extensions (patient route and disposable) with clamp, universal PVC connector and 1/4 adapter on the suction route.
5. Sterilized in ETO or available in a non-sterile version.
6. The bottle complies with ABNT NBR 12188:2016 standard (use between 11.81 inHg and 18.11 inHg) for vacuum pressure.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
15047	MULTIZAMM MAX 2000ml	S, V, Paciente 1/4"x 3,0m, Vácuo 1/4" x 3,0m, para lipo	10216350109	10
15052	MULTIZAMM MAX 3500ml	S, V, Paciente 1/4" x 2,50m, Vácuo 3/16" x 2,50m	10216350109	06
15152	MULTIZAMM MAX 3500ML NON STERILE	NS, V, Paciente 1/4" x 2,50m, Vácuo 3/16" x 2,50m	10216350110	06
15062	MULTIZAMM MAX 5000ML	S, V, Paciente 1/4" x 2,50m, Vácuo 3/16" x 2,50m	10216350109	05
15920	STERILE LID FOR MULTIZAMM PERDURA BOTTLE	S, V, Paciente 1/4" x 2,50m, Vácuo 3/16" x 2,50m	10216350118	30

*S Sterile*

*NS Non-sterile*

*V With Valve*

*Products subject to change without prior notice. Please contact us for availability.*