



## Indication

1. Extension for surgical procedures, gas conduction and aspiration.

## Description

1. Medium grade non-toxic transparent PVC extension cords in different diameters and lengths.
2. Flexible PVC connectors, graduated or smooth at the ends.
3. Latex-free.
4. Sterilized in ETO.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
9532	HOSPITAL EXTENSION CC 5,6 mm 2,0 m	S, 5,6mm, 2,0m, flexible connector	10216359002	100
9543	HOSPITAL EXTENSION CC 5,6 mm - 3,0 m	S, 5,6mm, 3,0m, flexible connector	10216359002	70
9544*	HOSPITAL EXTENSION CC 6,4 mm - 3,0 m	S, 6,4mm, 3,0m, flexible connector	10216359002	50
9553	HOSPITAL EXTENSION CC 5,6 mm - 6,0 m	S, 5,6mm, 6,0m, flexible connector	10216359002	35
9563	HOSPITAL EXTENSION CC 5,6 mm - 4,0 m	S, 5,6mm, 4,0m, flexible connector	10216359002	60
9753	HOSPITAL EXTENSION CC - 9,5 mm - 3,0 m - Lipo	S, 9,5mm, 3/8" 2,5mm x 3,0m, smooth connector	10216359002	20

*S Sterile*

*NS Non-sterile*

*V With Valve*

*Products subject to change without prior notice. Please contact us for availability.*