



Indication

1. PVC tube for connecting medical devices to oxygen lines.

Description

1. Non-toxic crystal green PVC extension, 3/16" in 1.5m or 2.0m lengths.
2. Latex-free.
3. Sterilized in ETO.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
9565	HOSPITALAR EXTENSION OX 4,7 mm - 1,5m	S, 4,7mm, 1,5m, smooth connector	10216359002	120
9566	HOSPITALAR EXTENSION OX 4,7 mm - 2,0m	S, 4,7mm, 2,0m, smooth connector	10216359002	110

S Sterile

NS Non-sterile

V With Valve

Products subject to change without prior notice. Please contact us for availability.