



PORTHOSP - REVZM08 | 02/2025

Indication

- 1.Perfusion of organs for transplantation.

Description

- 1.Two-way PVC syringe, 1/4" diameter, perforating tips plus protective caps and clamp, regulating flow cut-off in each of the ways.
- 2.Drip chamber and PVC connector at the end.
- 3.Latex-free.
- 4.Sterilized in ETO.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
8130	PERFUZAMM 2 WAY	S 1/4, extension 2,4m (approximately 2m + via serum 0.4)	10216350080	20

S Sterile
NS Non-sterile
V With Valve
Products subject to change without prior notice. Please contact us for availability.