



Indication

1. Pleural drainage with suction pressure regulation.
2. Helps control constant negative pressure.

Description

1. Made of rigid crystal PVC, guaranteeing resistance and transparency for safe monitoring.
2. 1,500ml capacity, with precise 50/50ml graduation and centimeter scale.
3. Equipped with PVC vacuum line, optimizing suction flow.
4. Contains PVC line with clamp and extension, making it easy to connect to collection bottles.
5. Sterilized in ETO.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
3415	PLEUROZAMM	E, SV, 3 vias, Paciente 1/4" x 2m, Vácuo 1/4" 1,0m	10216359004	10

S Sterile

NS Non-sterile

V With Valve

Products subject to change without prior notice. Please contact us for availability.