



## Indication

1. Aspiration of airway secretions or gastric/intestinal contents.
2. Use in procedures such as bronchoscopy, endoscopy and colonoscopy.

## Description

1. Rigid crystal PVC bottle, graduated every 100 ml, with a capacity of 1,000 ml.
2. Cap with internal thread and anti-reflux valve, equipped with bacteriological filter.
3. External thread compatible with vacuum gauge, ensuring a secure connection.
4. Crystal PVC extension, with clamp and protective cap.
5. Specific product, latex-free and ETO sterilized.
6. The bottle complies with ABNT NBR 12188:2016 standard (use between 11.81 inHg and 18.11 inHg) for vacuum pressure.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
1483	KIT STERZAMM 1000ml PLUS	S V, Pacient 3/16" x 1,0mm x 1,90m	10216350091	20

*\*Manometer for visualization. Not included.*

**S** Sterile

**NS** Non-sterile

**V** With Valve

*Products subject to change without prior notice. Please contact us for availability.*