



## Indication

1. Precise marking on patients' skin before surgical procedures.

## Description

1. Dye based on gentian violet, with antifungal and antibacterial properties, guaranteeing asepsis during the procedure.
2. Non-toxic and non-irritating, providing patient comfort and safety.
3. 1mm tip, ideal for fine, detailed strokes.
4. Comes with a 12 cm ruler, helping to ensure precise markings.
5. Sterilized by Gamma Ray, ensuring greater safety in use.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
10452	SURGICAL MARKING PEN	S, REGULAR, TIP1MM	-	25

S Sterile  
NS Non-sterile  
V With Valve  
Products subject to change without prior notice. Please contact us for availability.