



## Indication

1. Thoracic, mediastinal or pleural drainage.

## Description

1. Rigid crystal PVC bottle, graduated every 50ml, available in 500ml, 1,000ml and 2,000ml capacities.
2. Equipped with a quick-screw cap with 2 ports: one for the patient and one for the environment.
3. Patient route with 1.20m crystal PVC extension, choke, flow cut-off and plastic spiral to prevent elbowing.
4. Second route open to the environment, also functioning as a drain extension.
5. Sterilized in ETO.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
6200	TOROZAMM S 500ml PVC	S WATER SEAL, Pacient 5/16" x 1,8mm x 1,20m	10216350066	12
6400	TOROZAMM S 1000ml PVC	S, WATER SEAL, Pacient 5/16" x 1,8mm x 1,20m	10216350066	20
6601	TOROZAMM S 2000 ml PVC	S, WATER SEAL, Pacient 5/16" x 1,8mm x 1,20m	10216350066	10

*S Sterile*

*NS Non-sterile*

*V With Valve*

*Products subject to change without prior notice. Please contact us for availability.*