



Indication

- 1. Drainage and urine collection in a closed system, guaranteeing results in hourly diuresis control.
- 2. Helps monitor and measure bladder elimination, providing effective water control.

Description

- 1. Closed system hourly diuresis measurement set with graduated resin chamber (500ml) and collection bag (2000ml).
- 2. Integrated anti-reflux valve, which prevents urine from flowing back into the drip chamber, ensuring patient safety.
- 3. Flow controller for draining and emptying, positioned at the bottom of the drip chamber for better handling.
- 4. Urine sample collection device with calibrated connector, dispensing with the use of a needle for greater practicality.
- 5. Anti-traction clamp, which prevents the bladder catheter from being pulled by the weight of the system.
- 6. Latex-free and ETO sterilized.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
51100	UROZAMM 500 ADVANCED	S, rigid chamber 500ml, bag 2,000ml	10216350122	10

S Sterile
NS Non-sterile
V With Valve
Products subject to change without prior notice. Please contact us for availability.