



Indication

1. Drainage and urine collection in a closed system, guaranteeing results in hourly diuresis control.
2. Helps monitor and measure bladder elimination, providing effective water control.

Description

1. Closed system hourly diuresis measurement set with graduated resin chamber (500ml) and collection bag (2000ml).
2. Integrated anti-reflux valve, which prevents urine from flowing back into the drip chamber, ensuring patient safety.
3. Flow controller for draining and emptying, positioned at the bottom of the drip chamber for better handling.
4. Urine sample collection device with calibrated connector, dispensing with the use of a needle for greater practicality.
5. Anti-traction clamp, which prevents the bladder catheter from being pulled by the weight of the system.
6. Latex-free and ETO sterilized.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
5180	UROZAMM DH METER	S, rigid chamber 500ml, bag 2,000ml	10216350092	10

S Sterile

NS Non-sterile

V With Valve

Products subject to change without prior notice. Please contact us for availability.