



Indication

1. Drainage, collection and measurement in a closed system.

Description

- 1. Closed system urine collection bag, ensuring greater safety and hygiene.
- 2. 2,000ml capacity, with graduations starting at 25ml, allowing for precise accuracy.
- 3. Anti-reflux valve, which prevents urine from flowing back and reduces the risk of infection.
- 4. Protective cap, safe and hygienic comfort.
- 5. Latex-free and ETO sterilized.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
5195	UROZAMM STAR	S, BAG 2.000ml	10216359008	50

S Sterile
NS Non-sterile
V With Valve
Products subject to change without prior notice. Please contact us for availability.