



PORTHOSP - REVZM08 | 02/2025

Indication

- 1. Drainage, collection and measurement with barrier filter in closed system.

Description

- 1. Closed system urine collection bag, guaranteeing safety and control when monitoring diuresis.
- 2. 2,000ml capacity, graduated from 25ml, making it easy to measure the exact volume collected.
- 3. Pasteur chamber with integrated anti-reflux valve, reducing the risk of contamination and reflux.
- 4. T-shaped drain valve with cold guard, providing safe and hygienic emptying.
- 5. Latex-free and ETO sterilized.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
5190	UROZAMM TOP	S, BAG 2.000ml, Pasteur chamber, filter and waste valve with protector	10216350095	20

S Sterile  
NS Non-sterile  
V With Valve  
Products subject to change without prior notice. Please contact us for availability.