



PREMERE

Pressure Bag

The solution for fluids infusion

Adjustable diameter's bag for pressuring flexible bottles ou saline solutions bags.

- More resistant outer plastic film and thinner inner film allow faster and more consistent inflation.
- The bag surrounds the saline bottle 360°, allowing more even pressurization.
- 60mm diameter chrome-plate pressure gauge with easy-to-read scale.
- Velcro closure, positioned so as not to obstruct the view of the label, makes the installation of the serum bottle in the bag easier.
- Large inflation bulb and thick tube for faster inflation and deflation.
- Long term life time and easy to clean
- Available for 500ml or 1000ml bags



Code	Description	Anvisa	Box
7510	Bolsa Pressórica PREMERE 500ml	10216350119	1/6/16/48
7511	Bolsa Pressórica PREMERE 1000ml	10216350119	1/4/8/24



ZAMMI INSTRUMENTAL EIRELI

R. Bernardo Vasconcelos, 992 Parque Santa Lúcia, Duque de Caxias - RJ CEP: 25251-300

Tel: +55 21 2677-1500 | +55 11 5574-5767 | faleconosco@zammi.com.br

www.zammi.com.br