

Artrozamm

Indication

1. Artrozamm is a sterile product consisting of 1 set, with 2 or 4 extension lines, with varied lengths and diameters; drip chambers with protective caps, flow clamp on each line, and adapter for connection.
2. Product indicated for irrigation of cavities and/or perfusion of organs in surgeries, through the introduction of crystalloid solutions to clean the surgical cavity, to distend the cavity and keep it washed of blood and debris (particles), with controlled bleeding.

Descrição

1. PVC product with 2 or 4 extension ports, with diameters of 1/4".
2. Each port has a piercing tip, cap, and flow clamp.
3. Crystal PVC connector at the end, with option for luer lock connector
4. Latex-free product.
5. Sterilized in ETO.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
81003	ARTROZAMM 2 WAYS	S, 1/4" x 1,6mm x 2,0m	10216350113	25
81006	ARTROZAMM 4 WAYS	S, 1/4" x 1,6mm x 2,0m	10216350113	15
81007	ARTROZAMM 4 WAYS	S, 1/4" x 1,6mm x 2,0m	10216350113	15

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.

