


Avazamm

Indication

1. Suction of airway secretions through vacuum.

Description

1. Suitable system for secretion suction.
2. Available in 500 ml, 1000 ml and 2000 ml capacities.
3. Latex-free, minimizing allergy risks.
4. ETO sterilized or available in non-sterile version.
5. Integrated safety valve to prevent overflow.
6. The bottle complies with ABNT NBR 12188:2016.



CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
1425	KIT AVAZAMM 2000ml AR PLUS (KIT PARA LIPO)	NE, V, 1/4" x 1,6mm x 3,0m, 1/4" x 1,6mm x 3,0m	10216350070	10
1440	KIT AVAZAMM 1000ML AR PLUS	S, V, 3/16" x 1,0mm x 1,20m, 3/16" x 1,0mm x 0,50m	10216350091	20
1444	KIT AVAZAMM 1000ML AR PLUS	NE, V, 3/16" x 1,0mm x 1,20m, 3/16" x 1,0mm x 0,50m	10216350070	25
1492	KIT AVAZAMM 500ML AR PLUS	NE, V, 3/16" x 1,0mm x 1,20m, 3/16" x 1,0mm x 0,50m	10216350070	30
1495	KIT AVAZAMM 500ML AR PLUS	S, V, 3/16" x 1,0mm x 1,20m, 3/16" x 1,0mm x 0,50m	10216350091	25
1456	KIT AVAZAMM 1000ML AR PLUS	S, V, 3/16" x 1,0mm x 1,90m, 3/16" x 1,0mm x 0,50m	10216350091	20
1498	KIT AVAZAMM 500ML AR PLUS	S, V, 3/16" x 1,0mm x 1,90m, 3/16" x 1,0mm x 0,50m	10216350091	25

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.