

Drenosuczamm

Indication

1. Aspiration in large cavities, orthopedic surgeries, and organ retraction.
2. Conventional aspiration in the operating room.

Description

1. Rigid or flexible tip for suction, with flow regulator.
2. PVC extension of 3 meters, with flexible connector.
3. Non-toxic material.
4. Latex-free, ensuring greater safety.
5. Sterilized with ETO.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
1230	DRENOSUCZAMM EXTRA 20/8	S, 20/8 x 3,0m Rigid ti	10216350054	30
1222	DRENOSUCZAMM SOFT LONGO 18/7	S, 18/7 x 3,0m Rigid Tip	10216350054	30

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.