

# Extension for Surgical Center

## Indication

1. Extension for surgical procedures, gas conduction, and aspiration.

## Description

1. Transparent non-toxic crystal PVC extension, in different diameters and lengths.
2. Flexible PVC connectors, graduated or smooth at their ends.
3. Latex-free.
4. Sterilized in ETO.

CODE	DESCRIPTION	PRODUCT	ANVISA	BOX
9532	S, 5,6mm, 2,0m, flexible connect	Hospital Extension 5,6 mm 2,0 m	10216359002	100
9543	S, 5,6mm, 3,0m, flexible connect	Hospital Extension 5,6 mm - 3,0 m	10216359002	70
9544	S, 6,4mm, 3,0m, flexible connect	Hospital Extension 6,4 mm - 3,0 m	10216359002	50
9553	S, 5,6mm, 6,0m, flexible connect	Hospital Extension 5,6 mm - 6,0 m	10216359002	35
9563	S, 5,6mm, 4,0m, flexible connect	Hospital Extension 5,6 mm - 4,0 m	10216359002	60
9753	S, 9,5mm, 3/8" 2,5mm x 3,0m, smooth connect	Hospital Extension 9,5 mm - 3,0 m - Lipo	10216359002	20

*S Sterile*

*NE Non-sterile*

*V With Valve*

**Products subject to change without prior notice. Consult us about availability.**