

Extensão para Videolaparoscopia

Indication

1. Intermediate extension for the conduction of gases (CO₂; O₂...).

Description

1. Non-toxic clear PVC extension, with 1/4" and 3m in length.
2. Microbiological filter.
3. Latex-free.
4. Sterelized in ETO.



CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
8721	EXTENSÃO IVL	S, 1,4" x 1,6mm x 3,0m, FLEXIBLE CONNECTOR	10216350093	25

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.