





Irrigazamm Duo

Indication

1. Irrigation and suction in video surgery procedures.

Description

1. Double extension in non-toxic PVC of 1/4" x 3m.
2. Flexible PVC connector 1/4".
3. Flat choke for flow control.
4. Latex-free.
5. Sterilized in ETO.



CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
8200	IRRIGAZAMM DUO	S, 1/4" x 3,0m, WITH FLEXIBLE CONNECTOR, WITH 2 WAYS	10216359014	20

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.