

# Multizamm Max

## Indication

1. Aspiration of secretions and body fluids during surgical procedures.
2. Collection of secretions and purge through aspiration.

## Description

1. Rigid clear PVC bottle, graduated 100/100ml, with capacities of 2,000ml, 3,500ml, and 5,000ml.
2. Cap with quick screw with patient line identification and safety valve that stops suction when the bottle is full.
3. Anti-reflux valve with bacteriological filter, ensuring safety during use.
4. Two clear PVC extensions (patient line and vacuum) with flat clamp, flexible connector, and straight connector.
5. Sterilized with ETO.
6. The bottle complies with ABNT NBR 12188:2016 standard (use between 11.81 inHg and 18.11 inHg) of vacuum pressure.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
15047	MULTIZAMM MAX 2000ml	S, V, Pacient 1/4"x 3,0m, Vacuum 1/4" x 3,0m, para lipo	10216350109	10
15052	MULTIZAMM MAX 3500ml	S, V, Pacient 1/4" x 2,50m, Vacuum 3/16" x 2,50m	10216350109	06
15152	MULTIZAMM MAX 3500ML NAO ESTERIL	NE, V, Pacient 1/4" x 2,50m, Vacuum 3/16" x 2,50m	10216350110	06
15062	MULTIZAMM MAX 5000ML	S, V, Pacient 1/4" x 2,50m, Vacuum 3/16" x 2,50m	10216350109	05
15920	TAMPA ESTÉRIL PARA FRASCO MULTIZAMM PERDURA	S, V, Pacient 1/4" x 2,50m, Vacuum 3/16" x 2,50m	10216350118	30

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.