


Oxizamm

Indication

1. PVC tube to interconnect medical devices to oxygen lines.

Description

1. Transparent non-toxic PVC extension in green color, 3/16" in 1.5m or 2.0m length.
2. Latex-free.
3. Sterilized in ETO.



CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
9565	EXTENSÃO HOSPITALAR OX 4,7 mm - 1,5m	S, 4,7mm, 1,5m, connector smooth	10216359002	120
9566	EXTENSÃO HOSPITALAR OX 4,7 mm - 2,0m	S, 4,7mm, 2,0m, connector smooth	10216359002	110

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.