

Pleurozamm

Indication

1. Pleural drainage with suction pressure regulation.

Description

1. Manufactured in rigid crystal PVC.
2. 1,500 ml capacity with precise graduation.
3. Equipped with PVC vacuum line.
4. Contains PVC line with clamp and extension.
5. ETO sterilized and latex-free.



CODE	PRODUCT	DESCRIPTION	ANVISA	BOXCAIXA
3415	PLEUROZAMM	S, SV, 3 ways, Pacient 1/4" x 2m, Vácuo 1/4" 1,0m	10216359004	10

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.