



Sterzamm

Indication

1. Product intended for the collect of secretions by aspiration from the airways.
2. Use in procedures such as bronchoscopy, endoscopy, and colonoscopy.
3. The product meets the ABNT NBR 12188:2016 standard (use between 11.81 inHg and 18.11 inHg) for vacuum pressure.

Description

1. Rigid clear PVC bottle, graduated every 100 ml, with a capacity of 1,000 ml.
2. Cap with internal thread and anti-reflux valve, equipped with a bacteriological filter.
3. Rosca externa compatível com vacuômetro , garantindo conexão segura.
4. External thread compatible with vacuum gauge, ensuring a secure connection.
5. Specific product, latex-free and sterilized with ETO.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
1483	KIT STERZAMM 1000ml PLUS	E, V, Pacient 3/16" x 1,0mm x 1,90m	10216350091	20

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.