

Torozamm

Indication

1. Thoracic, mediastinal or pleural drainage.

Description

1. Rigid crystal PVC bottle available in 500 ml, 1,000 ml and 2,000 ml.
2. Internal tube with beveled terminal.
3. Patient line with crystal PVC extension.
4. ETO sterilized.



CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
6200	TOROZAMM S 500ml PVC	S, Watermark, Patient Copy 5/16" x 1,20m	10216350066	12
6400	TOROZAMM S 1000ml PVC	S, Watermark, Patient Copy 5/16" x 1,20m	10216350066	20
6601	TOROZAMM S 2000 ml PVC	S, Watermark, Patient Copy 5/16" x 1,20m	10216350066	10

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.