

Urozamm 500 Advanced

Indication

1. Product indicated for drainage and urine collection in a closed system, with hourly urine output monitoring.
2. There is no known contraindication for this product.

Description

1. Developed for urine collection and storage with maximum safety, ideal for hourly urine output monitoring.
2. It has 500ml in a rigid triple-graded chamber + 2000ml in a collection bag.
3. It features a foldable "T" tap and a bottom drain to facilitate emptying.



CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
51100	UROZAMM 500 Advanced	S, Rigid Chamber 500ml, bolsa 2.000ml	10216350122	10

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.