

Urozamm SF Star

Indication

1. Drainage, collection, and measurement in a closed system.

Description

1. Urine collector bag in a closed system, ensuring greater safety and hygiene.
2. Capacity of 2,000 ml, with graduation starting from 25 ml.
3. Anti-reflux valve, which prevents urine backflow and reduces the risk of infection.
4. Protective cap, safe and hygienic.
5. Latex-free and sterilized with ETO.



CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
5195	UROZAMM STAR	S, Bag 2.000ml	10216359008	50

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.