

Urozamm SF Top

Indication

1. The Urozamm SF-TOP is a product designed to collect and store urine from patients in a closed system, with a capacity of 2000 mL, graduated from 25 mL.
2. The product has a Pasteur chamber with an integrated anti-reflux valve and a "vent" for air entry with a bacteriological filter.

Description

1. Urine collector bag in a closed system, ensuring safety and control in monitoring urine output.
2. Capacity of 2,000ml, graduated from 25ml, facilitating accurate measurement of collected volume.
3. Pasteur chamber with integrated anti-reflux valve, reducing risks of contamination and reflux.
4. "T" drainage valve with cold protector, providing safe and hygienic emptying.
5. Latex-free and sterilized with ETO.

CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
5190	UROZAMM SF - TOP	Urine collector bag in a closed system S	10216350095	20

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.

