

Zammivac

Indication

1. Product indicated for continuous suction drainage.

Description

1. Postoperative vacuum suction aspirator.
2. Latex-free, offering greater safety for patients with sensitivities.
3. Sterilized with ETO.
4. Drains in diameters of 3.2mm, 4.8mm, and 6.4mm.
5. Includes universal Y connectors, clamps, and Zammivac connectors for easy adaptation and included.



CODE	PRODUCT	DESCRIPTION	ANVISA	BOX
1102I	ZAMMIVAC 3,2mm - 10FR - 1/8"	E, 3,2mm - 10FR - 1/8" NEEDLE 10FR, 20° E 28°	10216350056	12
1112I	ZAMMIVAC 4,8mm - 14FR - 3/16"	E, 4,8mm - 14FR - 3/16" NEEDLE 14FR, 20° E 28°	10216350056	12
1123I	ZAMMIVAC 6,4mm - 18FR - 1/4"	E, 6,4mm - 18FR - 1/4" NEEDLE 18FR, 20° E 28°	10216350056	12

S Sterile

NE Non-sterile

V With Valve

Products subject to change without prior notice. Consult us about availability.